

TiO₂ Anatase: What You Should Know But It's Not Always Explained

Text [Francesco Biaggini](#)

Anatase is one of the three main crystalline forms of titanium oxide (TiO₂), the other two being rutile and brookite. Among these, anatase is particularly valued for its unique photocatalytic and bacteriostatic properties, making it a material of great interest for biomedical and industrial applications. Anatase, characterized by a tetragonal structure, has demonstrated remarkable efficacy in the generation of reactive oxygen species (ROS) under UV irradiation, a mechanism that underlies its bacteriostatic capabilities.



In the field of biomedicine, anatase is mainly used in the coatings of medical devices, such as dental and orthopedic implants or high-end prosthetic abutments.

These coatings not only improve the biocompatibility of implants, but also impart bacteriostatic properties that help prevent post-operative infections, a significant problem in the clinical setting.

Anatase's ability to inhibit bacterial growth by reducing biofilm formation on implant surfaces makes it a material of choice for improving the safety and efficacy of medical devices.


One of the main approaches for the formation of anatase layers on metal surfaces is through titanium nitride (TiN) coating. Tin is known for its excellent hardness, wear resistance, and anti-corrosion properties, making it an ideal substrate for the deposition of anatase layers. Techniques such as physical vapor

deposition (PVD) and anodic oxidation are commonly used to create these coatings.

Anodizing is an electrochemical process that allows an oxide layer to form on the surface of titanium, which can later be converted into anatase. This method allows to precisely control the morphology and composition of the oxide layer, optimizing the photocatalytic and bacteriostatic properties of anatase.

In recent years, numerous studies have investigated the clinical applications of anatase, with a focus on dental and orthopedic implant coatings. Research continues to show that anatase coatings not only improve the biocompatibility and durability of implants, but also provide an effective barrier against bacterial infections. These developments represent a significant advance in the design of safer and longer-lasting medical devices.

In conclusion, anatase stands out as a highly valuable material due to its unique properties and applications in biomedical coatings. The combination of bacteriostatic properties and biocompatibility makes it a key component in improving the performance and safety of modern medical devices.



Bacteriostatic properties of anatase

Anatase is known for its exceptional bacteriostatic properties, which primarily stem from its ability to generate reactive oxygen species (ROS) when exposed to ultraviolet (UV) light. These reactive species include free radicals such as hydrogen peroxide (H₂O₂), hydroxyl radical (•OH), and singlet oxygen (¹O₂), which are highly reactive and can damage the cell membranes of bacteria, leading to their inactivation.

Mechanisms of action

Generation of Reactive Oxygen Species (ROS)

- Photocatalysis: Under UV irradiation, anatase catalyzes the formation of ROS through the excitation of electrons in the valence band of TiO₂, which jump into the conduction band, leaving behind positive gaps (holes). These electrons and holes

react with the water and oxygen present on the surface of the material, generating VSWR.

- Cell Membrane Damage: ROS attack the lipids of the bacterial cell membrane, causing lipid peroxidation, which leads to membrane disintegration and cell death.
- DNA damage: ROS can also penetrate inside bacterial cells and cause DNA damage, inhibiting cell replication and function.

Bacteriostatic effect

- Inhibition of Bacterial Proliferation: the action of ROS not only kills existing bacteria, but also inhibits their ability to proliferate, thus reducing the risk of biofilm formation. Bacterial biofilms are particularly difficult to treat in clinical settings because they protect bacteria from host immune defenses and antibiotics.
- Reduction of Antibiotic Resistance: The use of surfaces treated with anatase can reduce the need for antibiotics, thus decreasing the risk of bacterial resistance developing.

Scientific Evidence

Numerous studies have demonstrated the efficacy of anatase as a bacteriostatic agent. For example, in a study of anatase samples exposed to *Staphylococcus aureus* and *Escherichia coli*, a significant reduction in bacterial viability was observed compared to untreated samples. Samples treated with anatase showed up to a 99% reduction in bacterial growth under UV light conditions, confirming the efficacy of the ROS generated. Another study found that anatase coatings can prevent biofilm formation on medical surfaces, a common problem in dental and orthopedic implants. This effect is particularly important in the clinical setting, where infection prevention is crucial to the long-term success of implants.

In conclusion, anatase's bacteriostatic properties make it an ideal material for biomedical applications, especially in medical implant coatings. Its ability to generate ROS under UV irradiation provides an effective mechanism to inhibit bacterial growth and prevent infection, thereby helping to improve the safety and durability of medical devices.

Anatase Formation on Titanium Nitride Coatings

Titanium nitride (tin) is a material widely used as a coating to improve the surface hardness, wear resistance, and biocompatibility of medical devices. The addition of an anatase layer on tin can further enhance the bacteriostatic properties of the coating, combining the benefits of both materials. The formation of anatase on tin can be accomplished through several techniques, including physical vapor deposition (PVD) and anodic oxidation.

Deposition Processes

1. Physical Vapor Deposition (PVD)

- Working Principles: PVD is a technique that involves vaporizing the target material (in this case, titanium) in a reactive gas (nitrogen) environment to form a thin layer of titanium nitride on the surface of the substrate. Subsequently, an additional layer of anatase can be deposited by vaporization of TiO_2 in a controlled environment.
- Advantages: PVD offers high coating adhesion, good uniformity, and controllable composition. It is particularly useful for applications on complex surfaces.
- Limitations: The process requires expensive equipment and precise control of operating conditions to ensure coating quality.

2. Anodic oxidation

- Principles of Operation: Anodic oxidation is an electrochemical process that converts the surface of titanium into a layer of titanium oxide by the application of an electric current in a specific electrolyte. To form anatase, the anode process can be followed by heat treatment to achieve the desired crystalline phase.
- Critical Parameters: Key parameters include applied voltage, electrolyte composition, temperature, and treatment time. These factors influence the morphology, composition, and thickness of the oxide coating.
- Advantages: Anodic oxidation is relatively simple and can be performed at moderate temperatures. In addition, it allows precise control of the thickness of the oxide layer.
- Limitations: The quality of the coating can be affected by the purity of the substrate and the stability of the electrochemical process

Advantages of Titanium Nitride Coating

1. Corrosion Resistance

- TiN is highly resistant to corrosion, making it ideal for applications in wet or corrosive environments. This is especially important for dental and orthopedic implants, which must withstand chemical attack from body fluids.

2. Hardness and Wear Resistance

- TiN offers high surface hardness, which protects the substrate from scratches and wear. This extends the service life of the systems and reduces the need for frequent replacements.

3. Biocompatibility

- TiN is biocompatible and does not cause adverse reactions in human tissues, making it safe for use in medical devices.



Anatase, $\leq 10\text{nm}$


Anatase Formation on TiN

1. Deposition Techniques: After applying the tin coating, an anatase layer can be formed using techniques such as PVD or anodic oxidation followed by heat treatment.

The goal is to obtain a coating that combines the mechanical properties of tin with the bacteriostatic properties of anatase.

2. Parameter Optimization: The combination of deposition processes must be optimized to ensure good adhesion between layers and the formation of a crystalline anatase with high photocatalytic properties. This may include adjusting the temperature, treatment time, and electrolyte composition.

In summary, the formation of anatase coatings on titanium nitride substrates offers an advanced solution to improve the performance of medical implants. The combination of the properties of tin and anatase provides a durable, biocompatible surface capable of inhibiting bacterial growth, thus helping to reduce the risk of post-operative infections and improve the longevity of medical devices.



Anatase Formation on Oxidation

Anodizing is an electrochemical process that is widely used to improve the surface properties of metal materials, especially titanium.

This process not only improves the corrosion resistance and surface hardness of titanium, but also allows the formation of titanium oxides, including the anatase phase, which confers exceptional bacteriostatic properties. Anodizing is particularly effective in the production of highly controlled and uniform titanium oxide layers, which can subsequently be converted to anatase by specific heat treatment conditions.

Anodizing Process

1. Substrate Preparation

- Cleaning: The titanium substrate is thoroughly cleaned to remove any surface contaminants that may interfere with the anodizing process. Usually, acidic or alkaline solutions are used for this cleaning step.

- Activation: After cleaning, the substrate can be activated by chemical or heat treatments to improve the adhesion of the oxide coating.

2. Anodizing

- System Configuration: The titanium substrate is immersed in an appropriate electrolyte, such as phosphoric acid or sulfuric acid solutions, and connected as an anode in an electrical circuit.

- Application of Current: A constant or variable electric current is applied to the system, which induces the formation of an oxide layer on the surface of the titanium.

Anodizing parameters, such as voltage, current density, and electrolyte composition, are critical in determining the morphology and properties of the oxide film formed.

3. Anatase Formation

- Heat Treatment: The titanium oxide film formed during anodizing may contain amorphous phases or other crystalline forms of TiO_2 . To convert this film into anatase, a specific heat treatment is applied. The heating temperature, the treatment time and the environment (oxidizing or inert) play a crucial role in the formation of crystalline anatase.

Critical Parameters and Optimization

1. Voltage and Current Density

- Effects on Morphology: The applied voltage and current density directly affect the thickness and porosity of the oxide film. Higher voltages tend to produce thicker, more porous layers, while lower voltages result in thinner and more compact layers.

- Optimization for Anatase: To promote anatase formation, moderate voltages are often used that allow uniform growth of the oxide film.

2. Electrolyte Composition

- Types of Electrolytes: Commonly used electrolytes include phosphoric, sulfuric and boric acids, each with distinct effects on the composition and structure of the oxide film.

- Additives: Adding specific additives to the electrolyte can improve the quality of the oxide film, promoting anatase formation and improving bacteriostatic properties.

3. Heat Treatment

- Temperature and Time: Heat treatment for the conversion of anodic oxide to anatase generally takes place at temperatures between 300°C and 600°C. The holding time at these temperatures can vary from a few minutes to several hours, depending on the desired structure.

- Treatment Environment: The environment in which the heat treatment takes place (oxidizing or inert) can affect the purity and crystallinity of the anatase coating.

Practical Differences Between TiN and Anodizing from the Perspective of Durability, Stability and Effectiveness of the Anatase Layer

Mechanical Hardness

Titanium Nitride (TiN) coatings

- High Hardness: TiN is known for its exceptional hardness and wear resistance, which contribute to a long coating life. This makes it particularly suitable for applications that require high mechanical strength, such as dental and orthopedic implants.

- Corrosion Resistance: TiN offers excellent corrosion resistance, protecting the underlying metal substrate from corrosive body fluids, further extending the useful life of the medical device.

Anodization

- Variable Duration: The durability of anodized coatings depends on the quality of the oxide film formed and the exposure conditions. In general, anodized coatings are less resistant to mechanical wear than TiN, but they can offer good corrosion resistance.

- Maintenance: Anodized coatings may require periodic maintenance, especially in environments that are highly corrosive or subject to mechanical wear and tear.

Coating Stability over Time

Titanium Nitride (TiN) coatings

- Chemical Stability: TiN is extremely chemically stable, which means it does not degrade easily in the presence of harsh


chemicals. This makes it ideal for applications where long-term stability is critical.

- Strong Adhesion: TiN adheres very well to metal substrates, reducing the risk of delamination during prolonged use.

Anodization

- Modifiable Stability: The stability of anodized coatings can be improved through precise control of anodizing parameters and the application of post-anodizing treatments, such as heat treatment. However, they can be less chemically stable than TiN in extremely aggressive environments.

- Surface Morphology: The stability of the anodized coating can be affected by the porosity of the oxide film, which can make the coating more susceptible to mechanical and chemical damage if not properly treated.



Efficacy of Anatase Layer

Titanium Nitride (TiN) coatings

- Synergistic Effect: the application of an additional layer of anatase on TiN can significantly improve the bacteriostatic properties of the coating, combining the hardness and wear resistance of TiN with the photocatalytic properties of anatase.

- Bacteriostatic efficacy: When combined with anatase, TiN can offer excellent bacteriostatic efficacy, which is particularly useful for preventing infections in medical devices.

Anodization

- Photocatalytic Properties: Anodized coatings followed by heat treatment to form anatase exhibit excellent photocatalytic properties, making them highly effective in generating reactive oxygen species (ROS) that inhibit bacterial growth.

- Uniformity of the Coating: The anodizing technique allows the formation of uniform and well-adhering oxide layers, which can be optimized to improve bacteriostatic efficacy through the control of morphology and composition.

Final Confrontation

- Durability and Mechanical Strength: TiN offers superior durability and greater mechanical strength than anodized coatings. It is ideal for applications that require extreme robustness and a long service life.

- Chemical Stability: Both coatings are chemically stable, but TiN tends to be more resistant in aggressive environments.

- Bacteriostatic efficacy: While anodized coatings can offer excellent bacteriostatic properties due to the photocatalysis of anatase, the combination of TiN with anatase can provide an optimal solution that combines mechanical hardness and bacteriostatic activity.

In conclusion, the choice between TiN and anodizing depends on the specific application needs, balancing between durability, stability, and bacteriostatic properties to achieve the best possible result for the intended use of the medical device.

Clinical Applications and Recent Studies

Clinical Applications of Anatase in Implantable Medical Devices

The use of anatase as a coating for implantable medical devices has been shown to significantly improve the biocompatibility and bacteriostatic properties of these devices. Specifically, in dentistry, anatase finds application in dental implants, where infection prevention and promotion of bone integration are crucial for long-term success.



Dental Implants

- Infection Prevention: Anatase coatings on dental implants have been shown to significantly reduce the risk of post-operative infections. Anatase's ability to generate reactive oxygen species (ROS) under UV irradiation inhibits bacterial growth, reducing the formation of biofilms on the implant surface. This is especially important in the oral environment, where bacteria are abundant.

- Bone Integration: Anatase not only improves biocompatibility, but can also promote osseointegration, a critical process for successful dental implants. Studies have shown that anatase coatings can promote bone cell adhesion and proliferation, improving implant stability over time.

Orthopedic Prostheses

- Corrosion and wear resistance: In addition to dental implants, anatase is used in the coatings of orthopedic prostheses, such as hip and knee prostheses. These coatings not only improve the corrosion and wear resistance of the devices, but also provide a bacteriostatic barrier that reduces the risk of post-operative infections.

- Reduction of Post-Operative Complications: The bacteriostatic properties of anatase help prevent complications such as periprosthetic infections, one of the main causes of failure of orthopedic prostheses.



Recent Studies

1. Bacteriostatic Effect of Anatase-Coated Implants

- Studies reported on IJMS 2022: In a study presented in the paper "*Recent Advances in Metal-Based Antimicrobial Coatings for High-Touch Surfaces*", anatase coatings on titanium substrates were shown to exhibit significant antibacterial activity against various bacterial strains. The results showed a reduction in bacterial growth under conditions of exposure to UV light, confirming the effectiveness of the ROS generated.

2. Bone Integration and Cell Proliferation

- Research reported in the file "*Review of titanium surface modification techniques and coatings for antibacterial applications*": In this paper, anatase coatings have been shown to improve bone cell adhesion and proliferation. The study found that anatase-coated devices have better osseointegration than uncoated devices, improving the stability of dental and orthopedic implants in the long term.

3. Effects of anodizing on anatase

- Study reported in the file "*Titanium oxide antibacterial surfaces in biomedical devices*": This study explored the effect of anodizing on anatase formation and its bacteriostatic properties.

The results showed that anodized coatings with subsequent heat treatment to form anatase offer a significant reduction in bacterial growth and improved bone integration.

4. Anatase coatings on TiN

- Research reported in the file "*Antibacterial Applications of Anatase TiO₂ Nanoparticle*": This study investigated the combined use of tin coatings with anatase layers to improve the mechanical and bacteriostatic properties of medical devices. The results showed that the combination of Tin and anatase provides a highly durable and antibacterial surface, ideal for dental and orthopedic implant applications.

Conclusions and Discussions

The integration of anatase into implantable medical device coatings represents a significant step forward in improving the safety and efficacy of these devices. Its bacteriostatic and biocompatible properties not only help prevent infections, but also promote osseointegration, ensuring longer life and stability of implants. Recent studies support the use of anatase as an ideal coating material for a wide range of medical applications, particularly in dentistry and orthopedic surgery. But what does it really mean for the dentist or dental technician to consider the use of titanium nitride (tin) coatings or anodizing?

Clinical and Technical Considerations

1. Quality and Mechanical Characteristics

- TiN coatings are known for their exceptional hardness and wear resistance. This makes them particularly suitable for applications that require high mechanical strength, such as dental and orthopedic implants. The corrosion resistance of TiN protects the underlying metal substrate, ensuring a long life of the device.

2. Infection Prevention

- Anatase coatings, whether applied on tin or by anodizing, offer significant bacteriostatic properties due to the generation of reactive oxygen species (ROS) under UV irradiation. This bacteriostatic effect is crucial to reduce biofilm formation and prevent post-operative infections, which are particularly relevant in environments such as the oral one.

3. Biocompatibility and Osseointegration

- Anatase improves the biocompatibility of medical devices, promoting bone cell adhesion and proliferation. Studies show that anatase coatings can promote osseointegration, a critical process for the success of dental implants. This leads to greater stability and durability of the implant over time.

Stability Considerations for Anodized Coatings

Anodized coatings can offer good bacteriostatic properties, but they have some limitations compared to tin coatings:

- Anatase stability: Anatase layers formed by anodizing may be less stable than those on tins. In acidic environments, such as the oral cavity, thin layers of anodized oxide can dissolve more easily, compromising their long-term effectiveness.

- Lower Durability: The mechanical strength of anodized coatings is generally lower than that of TiN coatings. This can lead to faster wear and shorter device life, especially in high-stress applications such as dental implants.

- Impossible Maintenance: Maintenance of anodized coatings on implantable or dental devices is virtually impossible.

Anodized layers are typically very thin, a few microns, and are often used for aesthetic rather than functional purposes. Once the coating degrades, it cannot be easily restored.

Questions to Stimulate Reasoning

- How important is the durability and mechanical strength of the device to you?

While wear resistance and durability are priorities, TiN coatings offer a significant advantage due to their mechanical robustness.

- What is your main clinical concern: infection prevention or mechanical strength?

Both types of coating offer effective bacteriostatic properties, but anodized coatings can be less stable and long-lasting in acidic environments such as the mouth.

Stimulus to Discussion

- Is the technological compromise worth the cost savings?

Considering that anodized devices can be less expensive, it is important to consider whether the cost savings justify the potential trade-off in terms of mechanical strength and stability.

- How to assess the specific needs of the patient?

Every patient is unique, and needs can vary. For a patient with a high risk of infection, an anatase coating could offer a significant clinical benefit. For another patient, the durability and mechanical strength of the tin may be more crucial.

In conclusion, the choice of coating for implantable medical devices requires careful consideration of clinical needs, patient priorities, and specific operating conditions.

TiN coatings with anatase offer an ideal combination of mechanical strength and bacteriostatic properties, making them an excellent solution for many clinical applications. However, anodized coatings can be a viable alternative for specific applications, while still offering significant benefits in terms of infection prevention and biocompatibility, but with the awareness of their limitations in terms of stability in acidic environments and their inability to maintain.

As professionals, it is essential to carefully evaluate these factors to ensure the best possible outcome for patients.

Francesco Biaggini

